

Cutimed® Siltec® – The foam dressing that adapts to various exudate levels



Ordering Information

Cutimed® Siltec®	Pad size	Wound pad	Dressing size	Ref.-no.	Dressings/ box
Cutimed® Siltec® L (Light) Thin silicone foam dressing with Super Absorbers			5 cm x 6 cm	7328300	10
			10 cm x 10 cm	7328301	10
			15 cm x 15 cm	7328302	10
Cutimed® Siltec® Silicone foam dressing with Super Absorbers			5 cm x 6 cm	7328500	10
			10 cm x 10 cm	7328501	10
			10 cm x 20 cm	7328502	10
			15 cm x 15 cm	7328503	10
Cutimed® Siltec® B (Border) Self-adherent silicone bordered foam dressing with Super Absorbers	(4 x 4 cm) (6 x 6 cm) (8.2 x 8.2 cm) (6 x 16.5 cm)		7.5 cm x 7.5 cm	7328400	10
			10 cm x 10 cm	7328405	10
			12.5 cm x 12.5 cm	7328401	10
			10 cm x 22.5 cm	7328406	10
The anatomical specialists		Wound pad	Dressing size	Ref.-no.	Dressings/ box
Cutimed® Siltec® B Oval Self-adherent silicone bordered foam dressing with Super Absorbers			7 cm x 10 cm / oval	7328414	10
Cutimed® Siltec® Sacrum Self-adherent silicone bordered foam dressing for the sacral area with Super Absorbers			17.5 cm x 17.5 cm	7328700	5
			23 cm x 23 cm	7328701	5
Cutimed® Siltec® Heel Soft tack silicone foam dressing for the heel with Super Absorbers			16 cm x 24 cm / Heel 2D	7328805	5
			Heel 3D	7328601	5

Always follow the directions for use. These products are not available for purchase by the general public.

Cutimed® Siltec®


A Super Absorbent silicone foam dressing



Designed to work for you

The new Cutimed® Siltec® components

The right combination of specially designed materials makes **Cutimed Siltec** deliver the performance you expect from a foam dressing.



Cutimed® Siltec® range

Cutimed® Siltec®


When it comes to choosing a dressing that can be changed without damaging the wound or skin around it, **Cutimed® Siltec®** is driven by technology that meets the challenge. It is a soft and highly conformable foam dressing that adapts to all wound condition.

- Low to moderately exuding wounds
- Maintains critical moisture balance
- Skin protection with low adherence while reducing risk for maceration


Cutimed® Siltec® L

Cutimed® Siltec® L is a thinner version of **Cutimed® Siltec®** and is a soft and highly conformable foam dressing for low to moderate exuding wounds. **Cutimed® Siltec® L** can be changed without damage to the wound or surrounding skin.

- Soft and conformable wound dressing
- Protective dressing for compromised and fragile skin



Cutimed® Siltec®



Cutimed® Siltec® L

Cutimed® Siltec® B

Cutimed® Siltec® B has been designed to facilitate gentle adherence with effective exudate management for easy dressing application and removal.

Cutimed® Siltec® Sacrum

In order to ensure reliable exudate management when treating sacral ulcers, Essity has developed **Cutimed® Siltec® Sacrum**. This bordered foam dressing is particularly shaped to fit the awkward sacral area, offering an extra-thick wound pad for a high absorption capacity, plus providing additional cushioning.

- 5mm wound pad (with Super Absorbers) for moderate exudate levels and cushioning effect
- Wide, all-round adhering border
- Silicone wound contact layer for pain free and atraumatic removal
- Particularly shaped to fit the sacral region
- Shape allows for secure adherence with reduced friction and shear when immobile patient is repositioned
- Available in 2 sizes



Cutimed® Siltec® B



Cutimed® Siltec® Sacrum

Cutimed® Siltec® Heel

The dressing is designed to fit the anatomical shape of the heel and is folded around the heel at application.

- Super Absorbent soft silicone foam dressing
- Reduces pressure peaks, provides cushioning

Cutimed® Siltec® Heel 3D

Cutimed® Siltec® Heel 3D is especially designed to provide an optimal fit to every heel shape during ulcer management. Its anatomical shape provides cushioning and pressure relief from the heel to the Achilles tendon.

The extra thick foam of the dressing allows for reliable absorption of exudate.

- Super absorbent soft silicone foam dressing
- Ready-to-use, easy application



Cutimed® Siltec® Heel 2D



Cutimed® Siltec® Heel 3D